

For Research Use Only

# CoraLite® Plus 405 Anti-Mouse CD38 (90)



Catalog Number: **CL405-65059**

## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
CL405-65059	BC046312	Affinity purification
Size:	GeneID (NCBI):	CloneNo.:
100ug, 0.5 mg/ml	12494	90
Source:	UNIPROT ID:	Excitation/Emission maxima wavelengths:
Rat	P56528	399 nm / 422 nm
Isotype:	Full Name:	
IgG2a, kappa	CD38 antigen	

## Applications

Tested Applications:  
FC  
Species Specificity:  
Mouse

## Background Information

CD38, also known as ADP-ribosyl cyclase 1, is a type II transmembrane glycoprotein with a short N-terminal cytoplasmic tail, a single membrane-spanning domain, and a C-terminal extracellular region with four N-glycosylation sites (PMID: 2319135). The extracellular domain of CD38 has bifunctional enzyme activities that catalyze synthesis of cyclic ADP ribose from nicotinamide adenine dinucleotide (NAD) and hydrolysis of cyclic ADP ribose to adenosine diphosphoribose (PMID: 10636863). CD38 is expressed on a variety of hematopoietic and non-hematopoietic cells and is involved in diverse processes such as generation of calcium-mobilizing metabolites, cell activation, and chemotaxis (PMID: 25938500).

## Storage

Storage:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Anti-Mouse CD38 (CL405-65059, Clone:90) (red), or 0.5 ug Isotype Control. Cells were not fixed.