

For Research Use Only

CoraLite® Plus 405 Mouse IgG2a Isotype Control (MOPC-173)

Catalog Number: CL405-65244



Basic Information

Catalog Number:

CL405-65244

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Affinity purification

CloneNo.:

MOPC-173

Recommended Dilutions:

FC: 0.2 ug per 10⁶ cells in 100 µl suspension

Excitation/Emission maxima wavelengths:

399 nm / 422 nm

Source:

Mouse

Isotype:

IgG2a, kappa

Applications

Tested Applications:

FC

Species Specificity:

n/a

Positive Controls:

FC :

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 0.2 ug CoraLite® Plus 405 Anti-Mouse IgG2a Isotype Control (CL405-65244, Clone: MOPC-173) (red), CL405 Anti-Human CD5 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.