

For Research Use Only

CoraLite® Plus 405 Anti-Mouse CD206 Rabbit Recombinant Antibody

Catalog Number: CL405-98031



Basic Information

Catalog Number: CL405-98031	GenBank Accession Number: NM_008625.2	Purification Method: Protein A purification
Concentration: 100ug, 500 µg/ml	GeneID (NCBI): 17533	CloneNo.: 240344D12
Source: Rabbit	UNIPROT ID: Q2HZ94	Recommended Dilutions: FC: 0.25 ug per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: mannose receptor, C type 1	Excitation/Emission maxima wavelengths: 399 nm / 422 nm
Immunogen Catalog Number: EG0790	Calculated MW: 165 kDa	

Applications

Tested Applications: FC	Positive Controls: FC : mouse peritoneal macrophages,
Species Specificity: mouse	

Background Information

CD206, also named as MMR, CLEC13D and MRC1, is a type I membrane receptor that mediates the endocytosis of glycoproteins by macrophages. CD206 has been shown to bind high-mannose structures on the surface of potentially pathogenic viruses, bacteria, and fungi so that they can be neutralized by phagocytic engulfment. CD206 is a 170 kDa transmembrane protein which contains 5 domains: an amino-terminal cysteine-rich region, a fibronectin type II repeat, a series of eight tandem lectin-like carbohydrate recognition domains (responsible for the recognition of mannose and fucose), a transmembrane domain, and an intracellular carboxy-terminal tail. It is expressed on most tissue macrophages, in vitro derived dendritic cells, lymphatic and sinusoidal endothelia.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ FC receptor blocked mouse peritoneal macrophages were surface stained with CoraLite® Plus 488 Anti-Mouse CD107b / LAMP2 (ABL-93), and 0.25 ug CoraLite® Plus 405 Anti-Mouse CD206 Rabbit RecAb (CL405-98031, Clone:240344D12) or CoraLite® Plus 405 Rabbit IgG Isotype Control RecAb (CL405-98136, Clone: 240953C9). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



1x10⁶ FC receptor blocked mouse peritoneal macrophages were surface stained with 0.25 ug CoraLite® Plus 405 Anti-Mouse CD206 Rabbit RecAb (CL405-98031, Clone: 240344D12)(red) or CoraLite® Plus 405 Rabbit IgG Isotype Control RecAb (CL405-98136, Clone: 240953C9)(blue). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.