

For Research Use Only

CoraLite® Plus 488-conjugated TIMD4 Polyclonal antibody



Catalog Number: CL488-12008

Basic Information

Catalog Number:

CL488-12008

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2621

GenBank Accession Number:

BC008988

GeneID (NCBI):

91937

UNIPROT ID:

Q96H15

Full Name:

T-cell immunoglobulin and mucin
domain containing 4

Calculated MW:

378 aa, 42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human

Positive Controls:

IF : human tonsillitis tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488-conjugated TIMD4 antibody (CL488-12008) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488-conjugated TIMD4 antibody (CL488-12008) at dilution of 1:200.