

For Research Use Only



# CoraLite® Plus 488-conjugated ZWINT Polyclonal antibody

Catalog Number: **CL488-12282**

## Basic Information

Catalog Number:

CL488-12282

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2933

GenBank Accession Number:

BC020979

GeneID (NCBI):

11130

UNIPROT ID:

O95229

Full Name:

ZW10 interactor

Calculated MW:

277 aa, 31 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> HeLa cells were intracellularly stained with 0.8 ug CoraLite® Plus 488 Anti-Human ZWINT (CL488-12282) (red), or 0.8 ug Isotype Control. Cells were fixed and permeabilized with True-Nuclear Transcription Factor Buffer Set.