

For Research Use Only

CoraLite® Plus 488-conjugated Syntaxin 7 Polyclonal antibody



Catalog Number: **CL488-12322**

Featured Product

Basic Information

Catalog Number:

CL488-12322

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2975

GenBank Accession Number:

BC011975

GeneID (NCBI):

8417

UNIPROT ID:

O15400

Full Name:

syntaxin 7

Calculated MW:

261 aa, 30 kDa

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, canine

Positive Controls:

IF/ICC : A549 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol)
fixed A549 cells using CoraLite® Plus 488 Syntaxin
7 antibody (CL488-12322) at dilution of 1:200.