

For Research Use Only

CoraLite® Plus 488-conjugated MAGEA4 Polyclonal antibody

Catalog Number: CL488-12508



Basic Information

Catalog Number:

CL488-12508

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3193

GenBank Accession Number:

BC017723

GeneID (NCBI):

4103

UNIPROT ID:

P43358

Full Name:

melanoma antigen family A, 4

Calculated MW:

317 aa, 35 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ A431 cells were intracellularly stained with 0.4 ug CoraLite® Plus 488-conjugated MAGEA4 Polyclonal antibody (CL488-12508)(red), or 0.4 ug CoraLite® Plus 488-conjugated Rabbit IgG control Rabbit PolyAb (CL488-30000) (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).