

For Research Use Only

# CoraLite® Plus 488-conjugated MEK6 Polyclonal antibody



Catalog Number: CL488-12745

1 Publications

## Basic Information

Catalog Number:

CL488-12745

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3706

GenBank Accession Number:

BC012009

GeneID (NCBI):

5608

UNIPROT ID:

P52564

Full Name:

mitogen-activated protein kinase kinase 6

Calculated MW:

334 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC, FC (Intra)

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

IF/ICC : HepG2 cells,

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Hong Zhang	37690851	Adv Biol (Weinh)	IF

## Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488 MEK6 antibody (CL488-12745) at dilution of 1:200.



1X10<sup>6</sup> HepG2 cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human MEK6 (CL488-12745) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).