

For Research Use Only

# CoraLite® Plus 488-conjugated PDGFR beta Polyclonal antibody



Catalog Number: **CL488-13449**

Featured Product

## Basic Information

Catalog Number:

CL488-13449

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4285

GenBank Accession Number:

BC032224

GeneID (NCBI):

5159

UNIPROT ID:

P09619

Full Name:

platelet-derived growth factor receptor, beta polypeptide

Calculated MW:

1106 aa, 124 kDa

Observed MW:

160-180 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

