

For Research Use Only

CoraLite® Plus 488-conjugated PDGF-D Polyclonal antibody

Catalog Number: CL488-14075

Featured Product



Basic Information

Catalog Number:

CL488-14075

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5154

GenBank Accession Number:

BC030645

GeneID (NCBI):

80310

UNIPROT ID:

Q9GZP0

Full Name:

platelet derived growth factor D

Calculated MW:

43 kDa

Purification Method:

Antigen Affinity Purified

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, rat

Positive Controls:

IF/ICC : A549 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A549 cells using CoraLite® Plus 488 PDGFD antibody (CL488-14075) at dilution of 1:200.