

For Research Use Only

CoraLite® Plus 488-conjugated SMC3 Polyclonal antibody



Catalog Number: **CL488-14185**

Basic Information

Catalog Number:

CL488-14185

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5440

GenBank Accession Number:

BC047324

GeneID (NCBI):

9126

UNIPROT ID:

Q9UQE7

Full Name:

structural maintenance of chromosomes 3

Calculated MW:

142 kDa

Observed MW:

140 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

