

For Research Use Only

CoraLite® Plus 488-conjugated CRCP Polyclonal antibody



Catalog Number: **CL488-14348**

Basic Information

Catalog Number:

CL488-14348

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5689

GenBank Accession Number:

BC040107

GeneID (NCBI):

27297

UNIPROT ID:

O75575

Full Name:

CGRP receptor component

Calculated MW:

17 kDa

Observed MW:

20 kDa

Purification Method:

Antigen Affinity Purified

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Raji cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human CRCP (CL488-14348) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).