

For Research Use Only

CoraLite® Plus 488-conjugated PRMT3 Polyclonal antibody



Catalog Number: **CL488-17628**

Basic Information

Catalog Number:

CL488-17628

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11877

GenBank Accession Number:

BC037544

GeneID (NCBI):

10196

UNIPROT ID:

O60678

Full Name:

protein arginine methyltransferase 3

Calculated MW:

531 aa, 60 kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : HepG2 cells,

Background Information

PRMT3 (protein arginine methyltransferase 3) regulates protein functions by introducing asymmetric dimethylation marks at the arginine residues in proteins.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488 PRMT3 antibody (CL488-17628) at dilution of 1:200.