

For Research Use Only

# CoraLite® Plus 488-conjugated EDC4 Polyclonal antibody



Catalog Number: CL488-17737

Featured Product

## Basic Information

Catalog Number:

CL488-17737

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11784

GenBank Accession Number:

BC064567

GeneID (NCBI):

23644

UNIPROT ID:

Q6P2E9

Full Name:

enhancer of mRNA decapping 4

Calculated MW:

1401 aa, 152 kDa

Observed MW:

160 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:200-1:800

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, rat

Positive Controls:

IF : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 EDC4 antibody (CL488-17737) at dilution of 1:400, CL594-phalloidin (red).