

For Research Use Only

CoraLite® Plus 488-conjugated CHCHD7 Polyclonal antibody

Catalog Number: CL488-19911



Basic Information

Catalog Number:

CL488-19911

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13735

GenBank Accession Number:

BC002546

GeneID (NCBI):

79145

UNIPROT ID:

Q9BUK0

Full Name:

coiled-coil-helix-coiled-coil-helix
domain containing 7

Calculated MW:

97 aa, 12 kDa

Observed MW:

10 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

IF-P: 1:200-1:800

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF-P

Species Specificity:

human, mouse, rat

Positive Controls:

IF-P: mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse brain tissue using CoraLite® Plus 488 CHCHD7 antibody (CL488-19911) at dilution of 1:400. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).