

For Research Use Only

CoraLite® Plus 488-conjugated TAC1 Polyclonal antibody



Catalog Number: CL488-28599

Basic Information

| | | |
|---|--|---|
| Catalog Number: CL488-28599 | GenBank Accession Number: BC018047 | Purification Method: Antigen affinity purification |
| Size: 1000 µg/ml | GeneID (NCBI): 6863 | Recommended Dilutions: IF 1:50-1:500 |
| Source: Rabbit | UNIPROT ID: P20366 | Excitation/Emission maxima wavelengths: 493 nm / 522 nm |
| Isotype: IgG | Full Name: tachykinin, precursor 1 | |
| Immunogen Catalog Number: AG23477 | Calculated MW: 129 aa, 15 kDa | |

Applications

| | |
|--|--|
| Tested Applications: IF/ICC | Positive Controls: IF : PC-12 cells, |
| Species Specificity: Human, mouse, rat | |

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-12 cells using CoraLite® Plus 488 TAC1 antibody (CL488-28599) at dilution of 1:200.