

For Research Use Only

# CoraLite® Plus 488-conjugated IL-28A Monoclonal antibody



Catalog Number: CL488-60270

## Basic Information

<b>Catalog Number:</b> CL488-60270	<b>GenBank Accession Number:</b> BC113583	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 282616	<b>CloneNo.:</b> 3B4A1
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q8IZJO	<b>Recommended Dilutions:</b> IF 1:300-1:1200
<b>Isotype:</b> IgG1	<b>Full Name:</b> interleukin 28A (interferon, lambda 2)	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Immunogen Catalog Number:</b> AG18525	<b>Calculated MW:</b> 200 aa, 22 kDa	

## Applications

<b>Tested Applications:</b> IF-P	<b>Positive Controls:</b> IF : human tonsillitis tissue,
<b>Species Specificity:</b> human	

## Background Information

IL28A, also named as IFNL2 and ZCYT020, is a cytokine with immunomodulatory activity. It up-regulates MHC class I antigen expression. IL28A is a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IFNLR1. The ligand/receptor complex seems to signal through the Jak-STAT pathway. It plays a significant role in the antiviral immune defense in the intestinal epithelium.

## Storage

**Storage:**  
Store at -20°C. Avoid exposure to light.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 IL-28A antibody (CL488-60270, Clone: 3B4A1 ) at dilution of 1:600.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 IL-28A antibody (CL488-60270, Clone: 3B4A1 ) at dilution of 1:600.