

For Research Use Only

# CoraLite® Plus 488-conjugated CX3CL1 Monoclonal antibody



Catalog Number: **CL488-60339**

## Basic Information

Catalog Number:

**CL488-60339**

Size:

**1000 µg/ml**

Source:

**Mouse**

Isotype:

**IgG1**

Immunogen Catalog Number:

**AG0161**

GenBank Accession Number:

**BC001163**

GeneID (NCBI):

**6376**

UNIPROT ID:

**P78423**

Full Name:

**chemokine (C-X3-C motif) ligand 1**

Calculated MW:

**42 kDa**

Observed MW:

**90 kDa**

Purification Method:

**Protein G purification**

CloneNo.:

**2F5D5**

Excitation/Emission maxima wavelengths:

**493 nm / 522 nm**

## Applications

Tested Applications:

Species Specificity:

**human**

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

