

For Research Use Only

CoraLite® Plus 488-conjugated CX3CL1 Monoclonal antibody



Catalog Number: CL488-60339

Basic Information

Catalog Number:

CL488-60339

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0161

GenBank Accession Number:

BC001163

GeneID (NCBI):

6376

UNIPROT ID:

P78423

Full Name:

chemokine (C-X3-C motif) ligand 1

Calculated MW:

42 kDa

Observed MW:

90 kDa

Purification Method:

Protein G purification

CloneNo.:

2F5D5

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ PC-3 cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human CX3CL1 (CL488-60339, Clone:2F5D5) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).