

For Research Use Only

# CoraLite® Plus 488 Anti-Mouse CD81 (EAT2)



Catalog Number: CL488-65220

## Basic Information

Catalog Number:

CL488-65220

Size:

100ug, 500 µg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC011433

GeneID (NCBI):

12520

Full Name:

CD81 antigen

Purification Method:

Protein G purification

CloneNo.:

EAT2

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with 1 ug CoraLite® Plus 488 Anti-Mouse CD81 (CL488-65220, Clone:EAT2), or 1 ug Isotype Control. Cells were not fixed.