

For Research Use Only

# CoraLite® Plus 488 Anti-Rat CD43 (W3/13)



Catalog Number: **CL488-65237**

## Basic Information

Catalog Number:

CL488-65237

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

24796

UNIPROT ID:

P13838

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Affinity purification

CloneNo.:

W3/13

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

Rat

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 488 Anti-Rat CD43 (CL488-65237, Clone:W3/13), or 0.5 ug Control Antibody. Cells were not fixed.