

For Research Use Only

CoraLite® Plus 488-conjugated USP13 Monoclonal antibody

Catalog Number: CL488-66176

Featured Product



Basic Information

Catalog Number:

CL488-66176

Size:

1000 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10506

GenBank Accession Number:

BC016146

GeneID (NCBI):

8975

UNIPROT ID:

Q92995

Full Name:

ubiquitin specific peptidase 13
(isopeptidase T-3)

Calculated MW:

863 aa, 97 kDa

Observed MW:

97 kDa

Purification Method:

Protein G purification

CloneNo.:

6E8C10

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HepG2 cells,

Background Information

USP13(Ubiquitin carboxyl-terminal hydrolase 13) is also named as ISOT3 and belongs to the peptidase C19 family. It is involved in ubiquitin-dependent protein catabolism. USP13 is responsible for MITF deubiquitination, utilizing a short hairpin RNA library against known deubiquitinating enzymes. This protein has 2 isoforms produced by alternative splicing.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488-conjugated USP13 antibody (CL488-66176, Clone: 6E8C10) at dilution of 1:100.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 488 USP13 antibody (CL488-66176, Clone: 6E8C10) at dilution of 1:200.