

For Research Use Only

CoraLite® Plus 488-conjugated VPRBP Monoclonal antibody



Catalog Number: CL488-66392

Basic Information

| | | |
|--|--|---|
| Catalog Number: CL488-66392 | GenBank Accession Number: BC022792 | Purification Method: Protein A purification |
| Size: 1000 µg/ml | GeneID (NCBI): 9730 | CloneNo.: 1A7A8 |
| Source: Mouse | UNIPROT ID: Q9Y4B6 | Recommended Dilutions: IF 1:50-1:500 |
| Isotype: IgG2a | Full Name: Vpr (HIV-1) binding protein | Excitation/Emission maxima wavelengths: 493 nm / 522 nm |
| Immunogen Catalog Number: AG2184 | Calculated MW: 1506 aa, 169 kDa | |
| | Observed MW: 169 kDa | |

Applications

| | |
|---|--|
| Tested Applications: FC (Intra), IF/ICC | Positive Controls: IF : HepG2 cells, |
| Species Specificity: human | |

Background Information

VprBP was first identified as a protein that can interact with HIV-1 viral protein R (PMID: 11223251). It is a component of the CUL4A-RBX1-DDB1-VprBP/DCAF1 E3 ubiquitin-protein ligase complex that could interact with HIV-1 virus Vpr protein and HIV-2 virus Vpx protein (PMID: 18332868; 17314515; 18606781). VprBP is a 1,507-amino acid protein that contains conserved domains, including YXXY repeats, the Lis homology motif, and WD40 repeats. Through binding to Vpr, VprBP allows Vpr to modulate the catalytic activity of the CUL4-DDB1 complex, which in turn leads to the induction of G2 phase arrest in the virus-infected cells (PMID: 17630831). Recently it has been reported that VprBP is able to regulate the p53-induced transcription and apoptotic pathway (PMID: 22184063).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ K-562 cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human VPRBP (CL488-66392, Clone:1A7A8) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL488-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CoraLite® Plus 488 VPRBP antibody (CL488-66392, Clone: 1A7A8) at dilution of 1:200.