

For Research Use Only

CoraLite® Plus 488-conjugated TPRKB Monoclonal antibody



Catalog Number: CL488-66402

Basic Information

Catalog Number:

CL488-66402

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18770

GenBank Accession Number:

BC029492

GeneID (NCBI):

51002

UNIPROT ID:

Q9Y3C4

Full Name:

TP53RK binding protein

Calculated MW:

175 aa, 20 kDa

Observed MW:

18-20 kDa

Purification Method:

Protein G purification

CloneNo.:

2D5B9

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF : Jurkat cells,

Background Information

TPRKB, also named as CGI 121, belongs to the CGI121/TPRKB family. It interacts with TP53RK/PRPK and is widely expressed in all tissues.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using CoraLite® Plus 488 TPRKB antibody (CL488-66402, Clone: 2D5B9) at dilution of 1:200.