

For Research Use Only

CoraLite® Plus 488-conjugated PRDX3 Monoclonal antibody



Catalog Number: CL488-66810

Basic Information

Catalog Number:

CL488-66810

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG1062

GenBank Accession Number:

BC007062

GeneID (NCBI):

10935

UNIPROT ID:

P30048

Full Name:

peroxiredoxin 3

Calculated MW:

27 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

1F7F7

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ MCF-7 cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human PRX3 (CL488-66810, Clone:1F7F7) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL488-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).