

For Research Use Only

CoraLite® Plus 488-conjugated ADK Monoclonal antibody



Catalog Number: CL488-66929

Basic Information

Catalog Number: CL488-66929	GenBank Accession Number: BC003568	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 132	CloneNo.: 1G1A5
Source: Mouse	UNIPROT ID: P55263	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG1	Full Name: adenosine kinase	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG7466	Calculated MW: 41 kDa	

Applications

Tested Applications:
IF-P
Species Specificity:
Human, Mouse, Rat

Positive Controls:
IF : human liver cancer tissue,

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite® Plus 488 ADK antibody (CL488-66929, Clone: 1G1A5) at dilution of 1:200.