

For Research Use Only

# CoraLite® Plus 488-conjugated ADK Monoclonal antibody



Catalog Number: **CL488-66929**

## Basic Information

Catalog Number:

CL488-66929

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7466

GenBank Accession Number:

BC003568

GeneID (NCBI):

132

UNIPROT ID:

P55263

Full Name:

adenosine kinase

Calculated MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1G1A5

Recommended Dilutions:

IF-P 1:50-1:500

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF-P

Species Specificity:

Human, Mouse, Rat

Positive Controls:

IF-P : human liver cancer tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite® Plus 488 ADK antibody (CL488-66929, Clone: 1G1A5 ) at dilution of 1:200.