

For Research Use Only

CoraLite® Plus 488-conjugated TTPA Monoclonal antibody



Catalog Number: **CL488-67431**

Basic Information

Catalog Number:

CL488-67431

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG13587

GenBank Accession Number:

BC058000

GeneID (NCBI):

7274

UNIPROT ID:

P49638

Full Name:

tocopherol (alpha) transfer protein

Calculated MW:

32 kDa

Observed MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

3E7D10

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF-P

Species Specificity:

Human, Mouse, Rat, Pig

Positive Controls:

IF : human liver cancer tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite® Plus 488 TTPA antibody (CL488-67431, Clone: 3E7D10) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite® Plus 488 TTPA antibody (CL488-67431, Clone: 3E7D10) at dilution of 1:200.