

For Research Use Only

CoraLite® Plus 488-conjugated UBC12 Monoclonal antibody



Catalog Number:CL488-67482

Basic Information

Catalog Number:

CL488-67482

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5984

GenBank Accession Number:

BC058924

GeneID (NCBI):

9040

UNIPROT ID:

P61081

Full Name:

ubiquitin-conjugating enzyme E2M
(UBC12 homolog, yeast)

Calculated MW:

21 kDa

Observed MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1H8D2

Recommended Dilutions:

IF 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human, Mouse, Rat

Positive Controls:

IF : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using CoraLite® Plus 488 UBC12 antibody (CL488-67482, Clone: 1H8D2) at dilution of 1:200.