

For Research Use Only

CoraLite® Plus 488-conjugated ACF1 Monoclonal antibody



Catalog Number: **CL488-67947**

Basic Information

Catalog Number:

CL488-67947

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30607

GenBank Accession Number:

NM_013448

GeneID (NCBI):

11177

UNIPROT ID:

Q9NRL2

Full Name:

bromodomain adjacent to zinc finger domain, 1A

Calculated MW:

179 kDa

Purification Method:

Protein A purification

CloneNo.:

1C4F8

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-HUMAN ACF1 (CL488-67947, Clone:1C4F8) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).