

For Research Use Only

CoraLite® Plus 488-conjugated RAD23B Monoclonal antibody



Catalog Number: **CL488-67988**

Basic Information

Catalog Number: CL488-67988	GenBank Accession Number: BC020973	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 5887	CloneNo.: 2E6G9
Source: Mouse	UNIPROT ID: P54727	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG2a	Full Name: RAD23 homolog B (S. cerevisiae)	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG28658	Calculated MW: 409 aa, 43 kDa	
	Observed MW: 58 kDa	

Applications

Tested Applications: IF/ICC, FC (Intra)	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human, mouse, rat, pig	

Background Information

RAD23B, also named as p58, hHR23B and HR23B, belongs to the RAD23 family. It plays a central role both in proteasomal degradation of misfolded proteins and DNA repair. RAD23B is a central component of a complex required to couple deglycosylation and proteasome-mediated degradation of misfolded proteins in the endoplasmic reticulum that are retrotranslocated in the cytosol. Involved in DNA excision repair by stabilizing XPC protein. It may play a part in DNA damage recognition and/or in altering chromatin structure to allow access by damage-processing enzymes. RAD23B has two subunit: Heterodimer of a 125 kDa subunit (p125) and of a 58 kDa subunit (p58).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human RAD23B (CL488-67988, Clone:2E6G9) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL488-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).

Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 RAD23B antibody (CL488-67988, Clone: 2E6G9) at dilution of 1:200, CL594-Phalloidin (red).