

For Research Use Only

CoraLite® Plus 488-conjugated DHPS Monoclonal antibody



Catalog Number: CL488-67998

Basic Information

Catalog Number: CL488-67998	GenBank Accession Number: BC014016	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 1725	CloneNo.: 1E4F5
Source: Mouse	UNIPROT ID: P49366	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG1	Full Name: deoxyhypusine synthase	
Immunogen Catalog Number: AG31046	Calculated MW: 41 kDa	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
Human

Background Information

DHPS (Deoxyhypusine synthase) catalyzes the conversion of lysine to deoxyhypusine residue on the eukaryotic initiation factor 5A (eIF-5A) precursor using spermidine as the substrate. It belongs to the deoxyhypusine synthase family. This protein has 3 isoforms produced by alternative splicing.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug Coralite® Plus 488 Anti-Human DHPS (CL488-67998, Clone:1E4F5) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).