

For Research Use Only

CoraLite® Plus 488-conjugated SMS Monoclonal antibody



Catalog Number: **CL488-68040**

Basic Information

Catalog Number:

CL488-68040

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8840

GenBank Accession Number:

BC009898

GeneID (NCBI):

6611

UNIPROT ID:

P52788

Full Name:

spermine synthase

Calculated MW:

366 aa, 41 kDa

Observed MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1G8E9

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human

Positive Controls:

IF/ICC : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using Coralite® Plus 488 SMS antibody (CL488-68040, Clone: 1G8E9) at dilution of 1:100.