

For Research Use Only

CoraLite® Plus 488-conjugated SH3BP4 Monoclonal antibody



Catalog Number: **CL488-68045**

Basic Information

Catalog Number:

CL488-68045

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31238

GenBank Accession Number:

BC057396

GeneID (NCBI):

23677

UNIPROT ID:

Q9POV3

Full Name:

SH3-domain binding protein 4

Calculated MW:

963 aa, 108 kDa

Observed MW:

110 kDa

Purification Method:

Protein G purification

CloneNo.:

2G4H3

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

