

For Research Use Only

CoraLite® Plus 488-conjugated CDKN2A/P16-INK4A Recombinant monoclonal antibody

Catalog Number:CL488-81373-10



Basic Information

Catalog Number: CL488-81373-10	GenBank Accession Number: BC021998	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 1029	CloneNo.: 241438B2
Isotype: IgG	UNIPROT ID: P42771	Recommended Dilutions: IF/ICC: 1:50-1:500
Immunogen Catalog Number: AG1328	Full Name: cyclin-dependent kinase inhibitor 2A	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
	Calculated MW: 16 kDa	
	Observed MW: 16 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HEK-293 cells,
Species Specificity: human	

Background Information

P16-INK4A is also named as CDKN2A, MLM, Tumor suppressor ARF, Alternative reading frame. The tumor suppressor protein p16Ink4a (encoded from the CDKN2A locus) is often transcriptionally activated in cells undergoing senescence and is one of the main regulators of this program, and it is upregulated in multiple tissues during aging (PMID:17055429). p16-Ink4a is the principal member of the Ink4 family of CDK inhibitors. p16-Ink4a contributes to the regulation of cell cycle progression by inhibiting the S phase. p16Ink4a binds to CDK4/6, inhibiting cyclin D-CDK4/6 complex formation and CDK4/6-mediated phosphorylation of Rb family members. Expression of p16-Ink4a maintains the Rb family members in a hypophosphorylated state, which promotes binding to E2F1 and leads to G1 cell cycle arrest (PMID: 21297668).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using CoraLite® Plus 488 CDKN2A/P16-INK4A antibody (CL488-81373-10, Clone: 241438B2) at dilution of 1:200, CL594-phalloidin (red).