

For Research Use Only

CoraLite® Plus 488-conjugated BCL6 Recombinant antibody



Catalog Number: CL488-81754

Basic Information

Catalog Number:

CL488-81754

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15519

GenBank Accession Number:

BC150184

GeneID (NCBI):

604

UNIPROT ID:

P41182

Full Name:

B-cell CLL/lymphoma 6

Calculated MW:

706 aa, 79 kDa

Observed MW:

79-95 kDa

Purification Method:

Protein A purification

CloneNo.:

3K6

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human

Background Information

BCL6, a zinc finger transcription factor, contains an N-terminal BTB/POZ domain and C-terminal zinc finger DNA-binding motifs and represses transcription of a wide range of target proteins and microRNAs. BCL6 protein has been reported as a master regulator of B lymphocyte development and growth, and altered BCL6 protein expression was implicated in pathogenesis of diverse human hematologic malignancies, especially in the diffuse large B cell lymphoma (DLBCL). BCL6 is required for the development of T follicular helper T cells (TFH), a helper T cell subset required for the formation of mature and productive GCs. BCL6 has also been shown to play important regulatory roles in macrophages.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 Ramos cells were intracellularly stained with 0.8 μ g CoraLite® Plus 488 Anti-Human BCL6 (CL488-81754, Clone:3K6) (red), or 0.8 μ g Isotype Control. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).