

For Research Use Only

CoraLite® Plus 488-conjugated TBCD Recombinant antibody

Catalog Number: CL488-82911-7



Basic Information

Catalog Number:

CL488-82911-7

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6310

GenBank Accession Number:

BC003094

GeneID (NCBI):

6904

UNIPROT ID:

Q9BTW9

Full Name:

tubulin folding cofactor D

Calculated MW:

133 kDa

Observed MW:

130 kDa

Purification Method:

Protein A purification

CloneNo.:

230209F3

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : MCF-7 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed MCF-7 cells using CoraLite® Plus 488 TBCD antibody (CL488-82911-7, Clone: 230209F3) at dilution of 1:200.