

For Research Use Only

CoraLite® Plus 488-conjugated CD3 Delta Recombinant monoclonal antibody

Catalog Number: CL488-84214-4



Basic Information

Catalog Number: CL488-84214-4	GenBank Accession Number: BC070321	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 915	CloneNo.: 241523H2
Isotype: IgG	UNIPROT ID: P04234	Recommended Dilutions: IF-P: 1:50-1:500
Immunogen Catalog Number: EG2288	Full Name: CD3d molecule, delta (CD3-TCR complex)	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
	Calculated MW: 171 aa, 19 kDa	
	Observed MW: 20-25 kDa	

Applications

Tested Applications: IF-P	Positive Controls: IF-P: human tonsillitis tissue,
Species Specificity: human	

Background Information

CD3 is a complex of proteins that directly associates with the T cell receptor (TCR). The TCR/CD3 complex of T-lymphocytes consists of either a TCR alpha/beta or TCR gamma/delta heterodimer coexpressed at the cell surface with the invariant subunits of CD3 labeled gamma, delta, epsilon, zeta, and eta. The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction. CD3 is considered to be a pan-T cell marker for detection of normal and neoplastic T cells. CD3 delta/CD3D is a single-pass type I membrane protein which consists of an extracellular domain of 84 amino acids, a transmembrane region, and a cytoplasmic domain of 45 amino acids. Defects in CD3 delta are a cause of severe combined immunodeficiency.


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human tonsillitis tissue using CoraLite® Plus 488 CD3 Delta antibody (CL488-84214-4, Clone: 241523H2) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).