

For Research Use Only

CoraLite® Plus 488-conjugated FTSJ3 Recombinant antibody

Catalog Number: CL488-84311-6



Basic Information

Catalog Number:

CL488-84311-6

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9835

GenBank Accession Number:

BC036710

GeneID (NCBI):

117246

UNIPROT ID:

Q8IY81

Full Name:

FtsJ homolog 3 (E. coli)

Calculated MW:

847 aa, 97 kDa

Observed MW:

120-130 kDa

Purification Method:

Protein A purification

CloneNo.:

241507F2

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : MCF-7 cells,

Background Information

FTSJ3 is an RNA 2'-O-methyltransferase recruited by HIV to avoid innate immune sensing.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using Coralite® Plus 488 FTSJ3 antibody (CL488-84311-6, Clone: 241507F2) at dilution of 1:200.