

For Research Use Only

# CoraLite® Plus 488-conjugated CD44 Recombinant monoclonal antibody

Catalog Number:CL488-84369-5



## Basic Information

Catalog Number:

CL488-84369-5

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1458

GenBank Accession Number:

XM\_006234631.3

GeneID (NCBI):

25406

UNIPROT ID:

O70509

Full Name:

Cd44 molecule

Calculated MW:

40 kDa

Observed MW:

80 kDa

Purification Method:

Protein A purification

CloneNo.:

241648D1

Recommended Dilutions:

IF-P: 1:50-1:500

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF-P

Species Specificity:

rat

Positive Controls:

IF-P: rat spleen tissue,

## Background Information

CD44 is a type I transmembrane glycoprotein expressed on embryonic stem cells and in various levels on other cell types including connective tissues and bone marrow. CD44 expression is also upregulated in subpopulations of cancer cells and is recognized as a molecular marker for cancer stem cells (PMID: 29747682). It is a cell-surface receptor that mediates cell-cell and cell-matrix interactions through its affinity for hyaluronic acid (HA) and possibly also through its affinity for other ligands (PMID: 10694938). Adhesion with HA plays an important role in cell migration, tumor growth and progression. CD44 is also involved in lymphocyte activation, recirculation and homing, and in hematopoiesis.

## Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded rat spleen tissue using CoraLite® Plus 488 CD44 antibody (CL488-84369-5, Clone: 241648D1 ) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded rat spleen tissue using CoraLite® Plus 488 CD44 antibody (CL488-84369-5, Clone: 241648D1 ) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).