

For Research Use Only

CoraLite® Plus 488-conjugated THAP11 Recombinant antibody

Catalog Number: CL488-84373-6



Basic Information

Catalog Number:

CL488-84373-6

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19342

GenBank Accession Number:

BC012182

GeneID (NCBI):

57215

UNIPROT ID:

Q96EK4

Full Name:

THAP domain containing 11

Calculated MW:

314 aa, 34 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

241590E5

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using Coralite® Plus 488 THAP11 antibody (CL488-84373-6, Clone: 241590E5) at dilution of 1:200.