

For Research Use Only

CoraLite® Plus 488-conjugated IL-4RA/CD124 Recombinant antibody

Catalog Number: CL488-84414-6



Basic Information

Catalog Number: CL488-84414-6	GenBank Accession Number: NM_000418.4	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 3566	CloneNo.: 241824A11
Source: Rabbit	UNIPROT ID: P24394-1	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG	Full Name: interleukin 4 receptor	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: EG1684	Calculated MW: 90 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : Daudi cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Daudi cells using CoraLite® Plus 488 IL-4RA/CD124 antibody (CL488-84414-6, Clone: 241824A11) at dilution of 1:200.