

For Research Use Only

CoraLite® Plus 488-conjugated DOCK2-Specific Recombinant antibody

Catalog Number: CL488-84548



Basic Information

Catalog Number:

CL488-84548

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004946

GeneID (NCBI):

1794

UNIPROT ID:

Q92608

Full Name:

dedicator of cytokinesis 2

Calculated MW:

212 kDa

Purification Method:

Protein A purification

CloneNo.:

241350D8

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using CoraLite® Plus 488 DOCK2-Specific antibody (CL488-84548, Clone: 241350D8) at dilution of 1:200.