

For Research Use Only

CoraLite® Plus 488-conjugated TCB1 Recombinant monoclonal antibody

Catalog Number: CL488-84878-5



Basic Information

Catalog Number:

CL488-84878-5

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2437

GenBank Accession Number:

GeneID (NCBI):

UNIPROT ID:

P01852

Full Name:

Calculated MW:

19 kDa

Observed MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

241802B12

Recommended Dilutions:

IF-P: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF-P

Species Specificity:

mouse

Positive Controls:

IF-P : mouse thymus tissue,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse thymus tissue using CoraLite® Plus 488 TCB1 antibody (CL488-84878-5, Clone: 241802B12) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).