

For Research Use Only

CoraLite® Plus 488-conjugated Integrin Alpha V Recombinant monoclonal antibody



Catalog Number: **CL488-84883-5**

Basic Information

Catalog Number: CL488-84883-5	GenBank Accession Number: BC136442	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3685	CloneNo.: 242376D3
Isotype: IgG	UNIPROT ID: P06756	Recommended Dilutions: IF/ICC: 1:50-1:500
Immunogen Catalog Number: AG25825	Full Name: integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51)	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
	Calculated MW: 1048 aa, 116 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : A549 cells,
Species Specificity: human	

Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated alpha and beta subunits. Integrin alpha V, also known as CD51, is a type I transmembrane glycoprotein composed of an extracellular domain, a transmembrane region and cytoplasmic domain (PMID: 2430295; 1690718). It can form heterodimers with one of the five beta integrin subunits (beta 1, 3, 5, 6, 8) that recognize ligands containing the sequence Arg-Gly-Asp, such as vitronectin, fibronectin, and osteopontin (PMID: 37087799). Integrin alpha V mainly plays a role in extracellular matrix-mediated cell adhesion and migration, differentiation, wound healing, inflammation, tumorigenesis, proliferation, angiogenesis, and metastasis.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A549 cells using CoraLite® Plus 488 Integrin Alpha V antibody (CL488-84883-5, Clone: 242376D3) at dilution of 1:200, CL594-Phalloidin (red).