

For Research Use Only

CoraLite® Plus 488-conjugated TRIM27 Recombinant monoclonal antibody

Catalog Number: CL488-84927-6



Basic Information

Catalog Number:

CL488-84927-6

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2849

GenBank Accession Number:

BC013580

GeneID (NCBI):

5987

UNIPROT ID:

P14373

Full Name:

tripartite motif-containing 27

Calculated MW:

513 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

242538H5

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HepG2 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CoraLite® Plus 488 TRIM27 antibody (CL488-84927-6, Clone: 242538H5) at dilution of 1:200, CL594-Phalloidin (red).