

For Research Use Only

# CoraLite® Plus 488-conjugated SUGT1 Recombinant monoclonal antibody

Catalog Number:CL488-84963-5



## Basic Information

Catalog Number:

CL488-84963-5

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2275

GenBank Accession Number:

BC000911

GeneID (NCBI):

10910

UNIPROT ID:

Q9Y2ZO

Full Name:

SGT1, suppressor of G2 allele of SKP1  
(S. cerevisiae)

Calculated MW:

41 kDa

Observed MW:

38-41 kDa

Purification Method:

Protein A purification

CloneNo.:

242590G3

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : A431 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CoraLite® Plus 488 SUGT1 antibody (CL488-84963-5, Clone: 242590G3) at dilution of 1:200.