

For Research Use Only

CoraLite® Plus 488-conjugated TNPO1 Recombinant antibody

Catalog Number: CL488-85007-2

Featured Product



Basic Information

Catalog Number:

CL488-85007-2

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002270

GeneID (NCBI):

3842

UNIPROT ID:

Q92973

Full Name:

transportin 1

Calculated MW:

102 kDa

Purification Method:

Protein A purification

CloneNo.:

242269B11

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : U2OS cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite® Plus 488 TNPO1 antibody (CL488-85007-2, Clone: 242269B11) at dilution of 1:200, CL594-Phalloidin (red).



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using CoraLite® Plus 488 TNPO1 antibody (CL488-85007-2, Clone: 242269B11) at dilution of 1:200, CL594-Phalloidin (red).