

For Research Use Only

CoraLite® Plus 488-conjugated LGALS3BP Recombinant monoclonal antibody

Catalog Number: CL488-85883-5



Basic Information

Catalog Number: CL488-85883-5	GenBank Accession Number: NM_005567	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3959	CloneNo.: 250138C5
Isotype: IgG	UNIPROT ID: Q08380	Recommended Dilutions: IF/ICC: 1:50-1:500
Immunogen Catalog Number: EG3302	Full Name: lectin, galactoside-binding, soluble, 3 binding protein	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
	Calculated MW: 585 aa, 65 kDa	
	Observed MW: 97, 70 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HepG2 cells,
Species Specificity: human	

Background Information

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. Galectin-3-binding protein (LGALS3BP, also known as 90K or Mac-2 BP) is a secreted glycoprotein which binds galectin-3, galectin-1, beta1 integrins, collagens and fibronectin (PMID: 8034587; 8024581; 11146440; 9501082). It has been implicated in tumor metastatic processes, as well as in other cell adhesion and immune functions. Levels of LGALS3BP have been found elevated in the serum of patients with cancer and in those infected by the human immunodeficiency virus (HIV) (PMID: 7698018). Western analysis suggests that LGALS3BP is found in breast milk, semen, saliva, urine, and tears, in addition to serum (PMID: 8390986).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol)
fixed HepG2 cells using Coralite® Plus 488
LGALS3BP antibody (CL488-85883-5, Clone:
250138C5) at dilution of 1:200.