

For Research Use Only

CoraLite® Plus 488-conjugated RAP1GAP Recombinant monoclonal antibody

Catalog Number: **CL488-86011**



Basic Information

Catalog Number: CL488-86011	GenBank Accession Number: BC054490	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5909	CloneNo.: 25042288
Isotype: IgG	UNIPROT ID: P47736	Recommended Dilutions: IF/ICC: 1:50-1:500
Immunogen Catalog Number: AG5576	Full Name: RAP1 GTPase activating protein	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
	Calculated MW: 73 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : Jurkat cells, A549 cells
Species Specificity: human	

Background Information

Rap1 GTPase-activating protein (Rap1GAP) is an important tumor suppressor. Rap1GAP was associated with the AMPK and NF- κ B signaling pathways and that Rap1GAP inhibited AMPK/SIRT1 and activated the NF- κ B signaling pathway in model animals. Rap1GAP was verified to regulate the AMPK SIRT1/NF- κ B signaling pathway and exacerbate the damage to myocardial cells caused by ischemia and hypoxia (PMID: 38320624). Rap1GAP is involved in Ang II-induced cardiomyocyte hypertrophy through the regulation of autophagy (PMID: 33936386). Rap1GAP, which is the GTPase-activating protein of Rap1, can bind to Rap1, promote the transformation of GTP from its active to inactive form, and terminate signal transduction (PMID: 21295042).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using CoraLite® Plus 488 RAP1GAP antibody (CL488-86011, Clone: 250422B8) at dilution of 1:200.



Immunofluorescent analysis of (-20°C Methanol) fixed A549 cells using CoraLite® Plus 488 RAP1GAP antibody (CL488-86011, Clone: 250422B8) at dilution of 1:200.